

Supporting Information

G-Synth: Small Sequence Design and Extended Sequence Design for auditable synthesis-ready DNA and post-sequencing validation

Mohamed Merzoug^{a,*} · Zohra Yasmine Zater^{a, b} · Yasmine Saidi^a · Amaria Ilhem Hammadi^a · Marwa Aireche^a · Keltoum Bendida^a · Hadjer Soumia Bouderbala^a · Soheir Bouzidi^a · Djamal Saidi^a

^a Higher School of Biological Sciences of Oran, BP 1042 Saim Mohamed, Cité Emir Abdelkader, 31000 Oran, Algeria

^b Oran 1 University Ahmed Ben Bella, BP 1524 El M'Naouer, 31000 Oran, Algeria

* Corresponding author: Mohamed Merzoug; merzoug.mohamed@essb-oran.edu.dz

Contents: complete glargine sequences; SSD/ESD algorithmic invariants; reproducibility commands; unabridged experimental methods, synthesis script and QC record from the primary ACS draft; trace-file manifest; G-Synth evidence; secondary Geneious comparisons; and a scientific correction ledger.

S1. Exact glargine design inputs and outputs

Sequence type	Description	Sequence (preserved from primary draft)
Peptide Sequences	Insulin Glargine A-Chain	GIVEQCCTSICSLYQLENYCG
Peptide Sequences	Insulin Glargine B-Chain	FVNQHLCGSHLVEALYLVCGERGFFYTPKTRR
Optimized Nucleotide Sequences	Insulin Glargine A-Chain	GGCATCGTGGAACAGTGTGACACAGCATCTGCAG CCTGTACCAGCTGGAAAACACTGCGGCTAA
Optimized Nucleotide Sequences	Insulin Glargine B-Chain	TTTGTGAACAGCATCTGTGCGGCAGCCATCTGGT GGAAGCGCTGTACCTGGTGTGCGGCGAACGCGGC TTTTTTACACCCAAAAACCGCGCTAA

S1.1 Insulin glargine A-chain

Optimized coding input (66 bp):

5'→3' input

GGCATCGTGGAACAGTGTGACACAGCATCTGCAGCCTGTACCAGCTGGAAAACACTGCGGCTAA

Expected chain peptide: GIVEQCCTSICSLYQLENYCG

Forward order molecule (125 nt):

5'→3' forward

TATGGGTTCTTCTCACCACCACCACCACCTCTTCTGGTCTGGTGCCGCGTGGTTCTGGCATCGTGGAACAGTGTGCA
CCAGCATCTGCAGCCTGTACCAGCTGGAAAACACTGCGGCTAAC

Reverse order molecule (127 nt):

Validation flags: exact archived forward match; exact archived reverse match; exact duplex core; expected product translation.

S1.2 Insulin glargine B-chain

5'→3' input

TTTGTGAACCAGCATCTGTGCGGCAGCCATCTGGTGGAAGCGCTGTACCTGGTGTGCGGCGAACGCGGCTTTTTTTACA
CCCCAAAAACCCGCCGCTAA

5'→3' forward

TATGGGTTCTTCTCACCACCACCACCACCACTCTTCTGGTCTGGTGCCGCGTGGTTCTTTTGTGAACCAGCATCTGTGCGG
CAGCCATCTGGTGGAAGCGCTGTACCTGGTGTGCGGCGAACGCGGCTTTTTTACACCCAAAAACCCGCCGCTAAC

5'→3' reverse

TCGAGTTAGCGGCGGGTTTTGGGGTGAAAAAAGCCGCGTTCGCCGCACACCAGGTACAGCGCTTCCACCAGATGG
CTGCCGCACAGATGCTGGTTCACAAAAGAACCACGCGGCACCAGACCAGAAGAGTGGTGGTGGTGGTGAGAAGA
ACCCA

Validation flags: exact archived forward match; exact archived reverse match; exact duplex core; expected product translation.

Order molecule	Nucleotide sequence (5'→3')	Accession / translated product
Synthetic Insulin Glargine A-Chain Forward Sequence	TATGGGTTCTTCTACCAACCACCACCACTCTTCTGGTCT GGTGCCGCGTGGTTCTGGCATCGTGGAACAGTGTGCACC AGCATCTGCAGCCTGTACCAGCTGAAAACTACTGCGCT AAC	GenBank accession number PQ362993 MGSSHHHHHHSSGLVPRGSGIVEQCCTSICSLYQLENYCG (4352.87 Da)
Synthetic Insulin Glargine A-Chain Reverse Sequence	TCGAGTAGCCGAGTAGTITTTCCAGCTGGTACAGGCTGCA GATGCTGGTGACAGCACTGTTCCACGATGCCAGAACCAACG GGCACCAGACCAGAAGAGTGGTGGTGGTGGTGAGAA GAACCCA	GenBank accession number PQ362993 MGSSHHHHHHSSGLVPRGSGIVEQCCTSICSLYQLENYCG (4352.87 Da)
Synthetic Insulin Glargine B-Chain Forward Sequence	TATGGGTTCTTCTACCAACCACCACCACTCTTCTGGTCT GGTGCCGCGTGGTTCTTTTGGAACCAAGCATCTGCGCGC AGCCATCTGGTGGAAGCGCTGTACTGGTGTGCGGCGAAC GCGGCTTTTITACACCCCAAAAACCCGCGCTAAC	GenBank accession number PQ362994 MGSSHHHHHHSSGLVPRGSFVNQHLCSGSHLYEALYLCGER GFFYTPKTRR (5768.54 Da)

- Input sequences are normalized to unambiguous DNA and coding designs are translated before release.
- SSD releases one complementary forward/reverse synthesis pair; ESD tiles longer targets into multiple ordered complementary pairs with directional internal junctions.
- Recognition sequences and strand cut offsets determine terminal products; overhangs are not stored as pair-specific oligo literals.
- An enzyme supplies an initiating ATG only when the bases retained after cleavage place ATG at the coding boundary.
- The forward and reverse oligonucleotide cores must be exact reverse complements.
- Annealing and ordered ligation must reconstruct the designed molecule on both strands, base for base.
- Internal assembly overhangs must be unique against all used overhangs and their reverse complements; designs widen beyond four bases when required.
- Outer reconstructed ends must equal the selected enzyme products with correct sequence and polarity.
- A cloning vector must contain exactly one site for each selected enzyme; physical vector and insert ends must anneal.
- Coding-junction translation, stop-codon consequences, vector-tag outcomes and site regeneration are reported explicitly.
- Full sequence verification requires complete requested-region coverage and zero confident unexplained differences.

```
python tools/publication/design_glargine_case_study.py --manuscript <source-draft.docx> --output
publication_evidence/glargine_ab_design_and_cloning.json
python tools/publication/validate_insulin_correct_traces.py --root <reference-data-root> --trace-dir
publication_evidence/validated_traces --output-dir publication_evidence/sequencing_validation
python -m pytest gsynth_engine/tests -q
cd django_app && python -m pytest -q
cd ../frontend && npm test -- --run
```

Evidence object	Location	Purpose
Design/cloning JSON	publication_evidence/ glargine_ab_design_and_cloning.json	Exact molecules, segments, translations, junctions, warnings and source hash

Evidence object	Location	Purpose
Sequencing JSON	publication_evidence/ sequencing_validation/ glargine_approved_trace_validation.json	Approved-file policy, hashes, primary F/R consensus metrics, post hoc quality sensitivity and Biopython check
Design script	tools/publication/ design_glargine_case_study.py	Executable glargine case study with equality assertions
Trace script	tools/publication/ validate_insulin_correct_traces.py	Executable four-file sequencing analysis and figure generation

S4. Unabridged experimental oligonucleotide synthesis and cloning record

The record below was extracted from the primary ACS draft Manuscript_Merzoug_et_al_2025_ACS_--SB.docx (SHA-256 2338c8d46ecfd709ce4e1bf8007a752719fbd35795fb356b0dbfa953d187b810) without removing technical steps. Laboratory notebooks and instrument exports remain the authoritative operational record. Oligonucleotides were synthesized on a MerMade 4 instrument using 1- μ mol, 1000-Å universal supports and DMT-off processing.

Purification used an equal volume of cold butanol, 30 s vortexing, 30 min at -20°C and 15 min centrifugation at 14,000 rpm and 4°C . Pellets were washed twice with 1 mL 70% ethanol at -20°C , centrifuged for 10 min, air-dried and resuspended. Exact rotor radius is not recorded, so relative centrifugal force cannot be reconstructed from rpm; future protocols should report $\times g$.

Hybridization combined 10 μL of each complementary strand with 20 μL 2 $\times$ SSC containing 0.1% SDS. Samples were heated at 95°C for 5 min and cooled to room temperature over 30 min. pET-21a(+) digestion used 1 μg vector, 10 U NdeI, 10 U XhoI and 1 $\times$ CutSmart buffer in 50 μL for 2 h at 37°C . Gel-purified vector was ligated to duplex insert at a nominal 3:1 molar ratio with 50 ng vector and 400 U T4 DNA ligase in 20 μL at 16°C for 2 h followed by 4°C overnight. E. coli DH5 α cells were heat-shock transformed for 50 s at 42°C , recovered in SOC and selected on 100 $\mu\text{g}/\text{mL}$ ampicillin.

Measurement	A forward	A reverse	B forward	B reverse
ssDNA concentration (ng/ μL)	1702.52	1709.40	1688.51	1703.60
dsDNA concentration by ScanDrop (ng/ μL)	830	830	1936	1936
dsDNA A260/A280	1.92	1.92	1.88	1.88
Bioanalyzer apparent size (bp)	129	129	173	173
Bioanalyzer concentration (ng/ μL)	1.15	1.15	3.39	3.39

Caution: UV absorbance, gel mobility and Bioanalyzer sizing assess concentration, purity and approximate fragment size. They do not establish base identity. The different concentration values from ScanDrop and Bioanalyzer reflect method and dilution context and should not be combined without the original dilution records.




Figure S2. Primary wet-laboratory quality-control evidence. Agarose electrophoresis of the synthesized strands and hybridized duplexes, and Bioanalyzer electropherograms, are reproduced from the source experimental record. These assays support presence and approximate size, not nucleotide identity.

S5. Complete MerMade 4 synthesis program

Script detail	Reagent / manufacturer	No./V primes	Injection	Reaction	Equalize
INITIALIZATION	INITIALIZATION	INITIALIZATION	INITIALIZATION	INITIALIZATION	INITIALIZATION
3 × Washing	Anhydrous Acetonitrile, (emp BIOTECH, Berlin, Germany)	1 / 10	180 µL	00:00:02	00:00:12
RUN STEPS	RUN STEPS	RUN STEPS	RUN STEPS	RUN STEPS	RUN STEPS
Deblocking 1	3% Trichloroacetic Acid (TCA) in dichloromethane, (emp BIOTECH, Berlin, Germany)	1 / 10	130 µL	00:00:02	00:00:59
Deblocking 2	3% Trichloroacetic Acid (TCA) in dichloromethane, (emp BIOTECH, Berlin, Germany)	1 / 10	100 µL	00:00:02	00:00:45
2 × Washing	Anhydrous Acetonitrile, emp BIOTECH	1 / 10	180 µL	00:00:02	00:00:12
Coupling 1 (2 subinjections)	Activator reagent, (emp BIOTECH, Berlin, Germany)	1 / 10	70 µL	00:00:02	00:01:20
Coupling 1 (2 subinjections)	DNA Amidites (A, C, G, T), (Biosearch Technologies, Petaluma, CA, USA)	1 / 10	50 µL	00:00:02	00:01:20
Coupling 2 (2 subinjections)	Activator reagent, emp BIOTECH	1 / 10	50 µL	00:00:02	00:00:40
Coupling 2 (2 subinjections)	DNA Amidites (A, C, G, T), (Biosearch Technologies, Petaluma, CA, USA)	1 / 10	30 µL	00:00:02	00:00:40
Washing	Anhydrous Acetonitrile, (emp BIOTECH, Berlin, Germany)	1 / 10	180 µL	00:00:02	00:00:12
Capping (2 subinjections)	Cap A: Acetic anhydride in THF, emp BIOTECH	1 / 10	65 µL	00:00:02	00:00:50
Capping (2 subinjections)	Cap B: Acetic anhydride in THF, emp BIOTECH	1 / 10	65 µL	00:00:02	00:00:50

Script detail	Reagent / manufacturer	No./V primes	Injection	Reaction	Equalize
Washing	Anhydrous Acetonitrile, (emp BIOTECH, Berlin, Germany)	1 / 10	180 µL	00:00:02	00:00:12
Oxidation	Iodine in water/pyridine/THF (0.02 M), emp BIOTECH	1 / 10	80 µL	00:00:02	00:01:00
2 × Washing	Anhydrous Acetonitrile, (emp BIOTECH, Berlin, Germany)	1 / 10	180 µL	00:00:02	00:00:12
FINALIZATION	FINALIZATION	FINALIZATION	FINALIZATION	FINALIZATION	FINALIZATION
3 × Deblocking	3% Trichloroacetic Acid (TCA) in dichloromethane, (emp BIOTECH, Berlin, Germany)	1 / 10	150 µL	00:00:02	00:00:45
3 × Washing	Anhydrous Acetonitrile, (emp BIOTECH, Berlin, Germany)	1 / 10	180 µL	00:00:02	00:00:12
DEPROTECTION OUTSIDE THE OLIGOSYNTHESIZER	DEPROTECTION OUTSIDE THE OLIGOSYNTHESIZER	DEPROTECTION OUTSIDE THE OLIGOSYNTHESIZER	DEPROTECTION OUTSIDE THE OLIGOSYNTHESIZER	DEPROTECTION OUTSIDE THE OLIGOSYNTHESIZER	DEPROTECTION OUTSIDE THE OLIGOSYNTHESIZER
Cleavage & Deprotection	Concentrated ammonia solution, (Sigma-Aldrich, St. Louis, MO, USA)	/	1 mL	18 h at 70°C	/

S6. Primer record

Primer	Sequence (5'→3')	Reported T _m (°C)	Expected product (bp)
SpAB F	ATGGGTCTTCTCACCACCACCA	58.6	123 A / 156 B
SpA R	TTAGCCGAGTAGTTTTCCA	52.3	123 A
SpB R	TTAGCGCGGGTTTTTGG	53.8	156 B
SeqAB F	TGCATCCATATGGGTCTTCTCACCA CCACCA	67.2	144 A / 177 B
SeqA R	TGCATCCTCGAGTTAGCCGAGTAGT TTTCCA	65.8	144 A
SeqB R	TGCATCCTCGAGTTAGCGCGGGTT TTTGG	66.4	177 B
T7 F	TAATACGACTCACTATAGGG	52.0	307 A / 340 B
T7 R	GCTAGTTATTGCTCAGCGG	54.3	307 A / 340 B

The source table reported some A/B product sizes together in cells that do not apply to both reverse primers. The corrected product labels above associate SpA R/SeqA R with A and SpB R/SeqB R with B. Before a prospective run, each primer should be re-evaluated in G-Synth against the final recombinant plasmid and the laboratory's polymerase/buffer conditions.

S7. Wet-laboratory PCR and colony-screening evidence




Figure S3. PCR products obtained with insert-specific primers before cloning, reproduced from the primary experimental draft. SM, GeneRuler Express DNA Ladder; lanes 1-2, A insert; lanes 3-4, B insert.




Figure S4. Colony-screening PCR reproduced from the primary experimental draft. Lanes 1 and 3 are A and B products with insert-specific primers; lanes 2 and 4 are the corresponding products with T7 primers.

S8. Sequencing evidence manifest

The validation set is limited to A Forward Seq.ab1, A Reverse Seq.ab1, B Forward Seq.ab1 and B Reverse Seq.ab1. Their SHA-256 hashes are recorded in the evidence JSON; no other chromatogram contributes to any reported result. Each F/R pair was oriented and assembled before consensus coverage was evaluated. Consensus coverage, identity, bidirectional overlap and F/R agreement are the primary retrospective endpoints.

Chain	Ref. bp	Consensus cov.	Identity	F/R overlap	F/R agree	Consensus diff.	Biopython cov.
A	123	100.0%	100.0%	56.1%	100.0%	0	100.0%
B	156	100.0%	100.0%	71.2%	100.0%	0	100.0%




Figure S5. Primary G-Synth F/R consensus evidence map. The two oriented reads jointly cover each reference completely; green shading marks bidirectional overlap, with 100% agreement in both constructs.

S9. Supplementary interface evidence




Figure S6. Primer-template hybridization view. The annealing 3' segment is base-paired; the 5' cloning extension is explicitly unpaired in cycle 1 and appears in the predicted product after extension.




Figure S7. Detailed SSD duplex hybridization. The two order molecules form a 123-bp exact antiparallel core; the NdeI-derived 5'-TA and XhoI-derived 5'-TCGA terminal products remain exposed and are labelled as cohesive ends.

Using a pre-digested duplex, 125 bp with NdeI and XhoI ends. G-Synth will compare the observed strand geometry with the vector before enabling ligation.
Change insert or enzymes

3,490 bp 3,303 bp 123 bp 10 bp — aa

Compatibility checks

6 OF 6 PASSED

- ✓ **Overhangs are compatible**
NdeI and XhoI ends match the vector.
- ✓ **Both strands pair everywhere**
No mismatch in the insert duplex.
- ✓ **Each enzyme cuts the vector once**
Checked when the vector was cut; a second site is refused.
- ✓ **Orientation is forced**
The insert reads on the minus strand of the vector's own numbering, as every pET cassette does.
- ✓ **Sites are regenerated at both seams**
The insert can be cut back out — which is how the clone gets verified on a gel.
- ✓ **Reading frame survives the junction**
No reading frame was given.

Insert-vector end compatibility

VECTOR → DUPLEX INSERT → RECOMBINANT PRODUCT

VECTOR: **pET-21a(+)**
5,365 bp after digestion

INSERT: **construct**
125 bp duplex

PRODUCT: **Waiting for ligation**
review both junctions first

VECTOR → INSERT: **NdeI** (COMPLEMENTARY)
Vector end: 5'-TA-3' | | Insert end: 3'-AT-5'
5' cohesive ends align with opposite polarity.

INSERT → VECTOR: **XhoI** (COMPLEMENTARY)
Vector end: 5'-TCGA-3' | | | | Insert end: 3'-AGCT-5'
5' cohesive ends align with opposite polarity.

Every exposed insert end has the complementary vector end.
The button simulates joining the phosphodiester backbone; it does not represent an experimental ligation.

[Ligate compatible ends](#)

Notes

Figure S8. Pre-ligation vector-insert compatibility gate. The simple view preserves enzyme identity, overhang sequence and physical polarity while withholding product analyses.

Each enzyme cuts the vector once
Checked when the vector was cut; a second site is refused.

Orientation is forced
The insert reads on the minus strand of the vector's own numbering, as every pET cassette does.

Sites are regenerated at both seams
The insert can be cut back out — which is how the clone gets verified on a gel.

Reading frame survives the junction
No reading frame was given.

Insert-vector end compatibility

VECTOR → DUPLEX INSERT → RECOMBINANT PRODUCT

VECTOR: **pET-21a(+)**
5,365 bp after digestion

INSERT: **construct**
125 bp duplex

PRODUCT: **5,490 bp plasmid**
joined in silico

VECTOR → INSERT: **NdeI** (COMPLEMENTARY)
Vector end: 5'-TA-3' | | Insert end: 3'-AT-5'
5' cohesive ends align with opposite polarity.

INSERT → VECTOR: **XhoI** (COMPLEMENTARY)
Vector end: 5'-TCGA-3' | | | | Insert end: 3'-AGCT-5'
5' cohesive ends align with opposite polarity.

Every exposed insert end has the complementary vector end.
The button simulates joining the phosphodiester backbone; it does not represent an experimental ligation.

[Ends ligated](#)

Map

[Circular](#) [Linear](#) [Both](#)

☒ Restriction sites ☐ Only the pair used for cloning ☐ Include multi-cutters (all catalogued sites)

17 SITES DRAWN · SINGLE-CUTTERS · CLONING PAIR · ORIGIN-CROSSING SITE SPLIT ACROSS BOTH ENDS

Map showing restriction sites: XhoI, 6xHis, T7 terminator, BgIII, +7.

Figure S9. Product state after explicit in-silico ligation. The recombinant length and downstream map are displayed only after both junctions have passed.




Figure S10. HindIII-containing diagnostic digest and predicted gel. The restriction site and complete digest fragments are displayed with a named broad-range marker. The panel is labelled as an in-silico prediction.




Figure S11. Restriction map showing HindIII as a named site derived from cut geometry.

S10. Secondary Geneious comparison and independent confirmation

The supplied experimental draft used Benchling for reverse complements and Geneious Prime 2024.0.7 for cloning views and trace alignments. G-Synth independently regenerated the same four synthesis molecules from the two optimized coding inputs, and its recombinant maps placed the inserts downstream of the pET-21a(+) T7/lac expression elements with preserved coding frames. The archived Geneious figures and G-Synth therefore agree qualitatively on construct identity and placement. A new automated Geneious run was not performed

because the reproducibility environment did not include a licensed command-line execution path. Biopython 1.87 was used as the independent, executable alignment comparator.

A prospective head-to-head study should freeze one set of reference plasmids and raw trace files, define identical trimming thresholds and variant rules, blind sample labels, and compare: read placement; orientation calls; covered intervals; mismatches/indels; consensus; runtime; manual interventions; exportability; and final verdict. The present study is a functional cross-check, not a statistical superiority benchmark.




Figure S12. Archived Geneious Prime 2024.0.7 duplex representation of NdeI/XhoI-compatible A- and B-chain inserts. This is secondary qualitative concordance; G-Synth SSD is the primary in-silico reconstruction.




Figure S13. Archived Geneious Prime construct annotation for the glargine A- and B-chain pET-21a(+) designs. G-Synth generated the primary coordinate, ORF and junction validation reported in the main article.




Figure S14. Archived Geneious Prime alignment made with the same four author-approved construct-confirmation files used by G-Synth. It provides secondary visual concordance and was not treated as an independent blinded benchmark.

S11. Scientific correction ledger for the supplied draft

Source-draft implication	Correction in the G-Synth manuscript	Reason
Sequencing showed complete 100% identity	Oriented F/R assembly produced 100% consensus coverage and identity; overlap support was 56.1% (A) and 71.2% (B), with 100% F/R agreement	Separates assembled coverage from bidirectional confirmation without obscuring the complete consensus
The method eliminates enzymatic digestion	Insert preparation is PCR-free and pre-engineered; the pET-21a(+) vector is still digested with NdeI/XhoI	Vector digestion is experimentally required
The constructs enable insulin glargine production	Construct design and cloning were validated; expression, cleavage, refolding and bioactivity were not tested	DNA construction does not establish functional insulin
His-tags and linkers ensure folding remains unaffected	Sequence presence and tag position are reported; folding and cleavage remain experimental	No protein-level evidence was supplied

Source-draft implication	Correction in the G-Synth manuscript	Reason
The workflow is broadly scalable and cost-effective	Potential advantages are framed as hypotheses for prospective benchmarking	No controlled cost, yield or scale comparison was performed
Geneious is the primary validator	G-Synth is the primary design/cloning/trace system; archived Geneious views and executable Biopython provide independent confirmation	Matches the software paper's validation hierarchy
One oligo length per insert	Forward/reverse lengths are reported separately (125/127 nt and 158/160 nt)	Cohesive-end geometry is strand-asymmetric

S12. Release and submission checklist

- Obtain written approval from every listed author for author order, affiliations, CRediT roles and AI disclosure.
- Run prospective bidirectional sequencing until both junctions and 100% of each insert are covered at the predefined quality threshold.
- Deposit raw trace files, reference FASTA/GenBank files, evidence JSON and analysis scripts in a stable repository, subject to institutional approval.
- Publish the version 1.0.0 GitHub release, connect Zenodo and insert the versioned DOI in the manuscript and CITATION.cff.
- Perform a clean-clone execution of all three test suites and both publication scripts; retain logs and environment versions.
- Prepare the target journal's graphical summary, high-resolution figures and machine-readable constructs; confirm its current author guidelines and retain the limitation that functional insulin glargine was not demonstrated.